

NEW FROM TETRACORE®: T-COR 4™ REAL-TIME PCR THERMOCYCLER



Thermocycler Highlights:

- Highly portable – four sample wells
- Multiplex capability – 2 color + IC, or 3 color
- Rechargeable battery
- Seven hour thermocycling battery life
- Runs in stand-alone mode or PC daisy chain

Tetracore® Real-Time PCR reagents and portable PCR thermocycler analyze unknown samples for the rapid confirmatory analysis of potential BW agents. Ideal for use by law enforcement, local, state, and federal governments.



We manufacture dried-down (room temperature) real-time PCR reagents for:

<i>B. anthracis</i> Lethal Factor (pXO1)	Orthopox Virus	Burkholderia	Classical Swine Fever Virus
<i>B. anthracis</i> capA (pXO2)	Ricin	Rinderpest	African Swine Fever Virus
<i>Brucella</i> spp.	Botulinum toxin	Capripox	Venezuelan Equine Encephalitis Virus
<i>F. tularensis</i>	SEB	Peste des Petits Ruminants	Foot and Mouth Disease Virus
<i>Y. pestis</i>			Japanese Encephalitis Virus



NEW FROM TETRACORE®: T-COR 4™ REAL-TIME PCR THERMOCYCLER



Thermocycler Highlights:

- Highly portable – four independent sample wells
- Multiplex capability – 1 color + IC, or 2 color
- Rechargeable battery
- Eight hour thermocycling battery life
- Runs in stand-alone mode or PC daisy chain

Tetracore® Real-Time PCR reagents and portable PCR thermocycler analyze potential threats for the rapid confirmatory analysis of potential BW agents and high-impact veterinary diseases.



We manufacture dry (room temperature) real-time PCR reagents* for:

B. anthracis Lethal Factor (pXO1)
B. anthracis capA (pXO2)
Brucella spp.
Y. pestis
F. tularensis

Orthopox Virus
Botulinum toxin
Ricin
SEB

*for research & development use only

Burkholderia
Rinderpest
Capripox
Peste des Petits Ruminants

Classical Swine Fever Virus (USDA Lic.)
African Swine Fever Virus
Foot and Mouth Disease Virus
Japanese Encephalitis Virus
Venezuelan Equine Encephalitis Virus



Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

USDA Licensed VetAlert Real Time PCR Kits				
Product Name	Catalog Number	Product Description	Test Format & Instrument Platforms	Unit Size
CSF RT-PCR (SC)	TC-9002-064	Diagnosis for Classical Swine Fever (CSF) Disease in Swine	Dry Reagents for use with: Cepheid SmartCycler®	1 Kit (64 reactions)
CSF RT-PCR (LC)	TC-9019-064	Diagnosis for Classical Swine Fever (CSF) Disease in Swine	Dry Reagents for use with: Roche LightCycler®, ABI Prism™, BioRad iCycler®, Stratagene®	1 Kit (64 reactions)
Johne's Real Time PCR *	TC-9828-100	Diagnosis for Johne's Disease in Cattle	Wet Reagents for use with: Cepheid SmartCycler®, Roche Light Cycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	1 Kit (100 reactions)
Johne's Combo-Pak*	TC-9011-100	Diagnosis for Johne's Disease in Cattle plus fecal DNA extraction	Wet Reagents for use with: Cepheid SmartCycler®, Roche Light Cycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	1 Kit (100 reactions + 100 extractions)
VetAlert Extraction System				
Product Name	Catalog Number	Product Description	Test Format & Instrument Platforms	Unit Size
MAP Extraction System*	TC-9014-100	Kit For the Extraction of Fecal DNA	N/A	1 Kit (100 reactions)

***Contact Tracy Fecteau at 240-268-5400 or tfecteau@tetracore.com**

Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

Veterinary Real Time PCR Reagents				
Product Name	Catalog Number	Product Description	Test Format & Instrument Platforms	Unit Size
NextGen PRRSV Multiplex RT-PCR	TC-9039-096 TC-9039-192	Reagents for the simultaneous <u>detection & differentiation</u> of both North American and European-like strains of Porcine Reproductive and Respiratory Syndrome (PRRS) Viral RNA	Wet reagents for use with: Cepheid SmartCycler®, Roche LightCycler®, ABI PRISM™, Stratagene®	6 Vials (96 reactions for TC-9039-080) 4 Vials (192 reactions for TC-9039-240)
EURO PRRSV Multiclade RT-PCR	TC-9034-096	Reagents for the specific detection of European-like/Lelystad Strains of Porcine Reproductive and Respiratory Syndrome (PRRS) Viral RNA	Wet reagents for use with: Cepheid SmartCycler®, Roche LightCycler®, ABI PRISM™, Stratagene®	6 Vials (96 reactions)
NA PRRSV RT-PCR	TC-9013-096	Reagents for the specific detection of North American Strains of Porcine Reproductive and Respiratory Syndrome (PRRS) Viral RNA	Wet reagents for use with: Cepheid SmartCycler®	6 Vials (96 reactions)
FMD RT-PCR (LC)	TC-9015-064	Reagents for specific detection of Foot and Mouth Disease (FMD) Viral RNA	Dry reagents for use with: Roche LightCycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	1 Unit (64 reactions)
FMD RT-PCR (SC)	TC-9006-064	Reagents for specific detection of Foot and Mouth Disease (FMD) Viral RNA	Dry reagents for use with: Cepheid SmartCycler®	1 Unit (64 reactions)
WNV RT-PCR	TC-9003-080	Reagents for specific detection of West Nile Virus (WNV) RNA	Wet reagents for use with: Cepheid SmartCycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	5 Vials (80 reactions)
ASF PCR (LC)	TC-9017-064	Reagents for the specific detection of African Swine Fever (ASF) Viral DNA	Dry reagents for use with: Roche LightCycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	32 vials (64 reactions)

9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
 Phone: 240.268.5400 * Fax: 240.268.1107
 ® Tetracore is a Registered Trademark

Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

Veterinary Real Time PCR Reagents				
Product Name	Catalog Number	Product Description	Test Format & Instrument Platforms	Unit Size
ASF PCR (SC)	TC-9038-064	Reagents for the specific detection of African Swine Fever (ASF) Viral DNA	Dry reagents for use with: Cepheid SmartCycler®	1 Unit (64 reactions)
Capripox Virus PCR (LC)	TC-9029-064	Reagents for the specific detection of Capripox Viral DNA	Dry reagents for use with: Roche LightCycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	32 vials (64 reactions)
Rinderpest Virus RT-PCR (LC)	TC-9030-064	Reagents for the specific detection of Rinderpest Viral RNA	Dry reagents for use with: Roche LightCycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	32 vials (64 reactions)
Bluetongue Virus RT-PCR (LC)	TC-9031-064	Reagents for the specific detection of Bluetongue Viral RNA	Dry reagents for use with: Roche LightCycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	32 vials (64 reactions)
Peste des Petits Ruminants Virus RT-PCR (LC)	TC-9032-080	Reagents for the specific detection of Peste des Petits Ruminants (PPR) Viral RNA	Dry reagents for use with: Roche LightCycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	5 Vials (80 reactions)
Burkholderia PCR (LC)	TC-9033-064	Reagents for the specific detection of Burkholderia DNA	Wet or Dry reagents for use with: Roche LightCycler®, ABI PRISM™, BioRad iCycler®, Stratagene®	32 vials (64 reactions)
Burkholderia PCR (SC)	TC-9041-064	Reagents for the specific detection of Burkholderia DNA. Includes Inhibition Control.	Dry reagents for use with: Cepheid SmartCycler®	1 Unit (64 reactions)

9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
 Phone: 240.268.5400 * Fax: 240.268.1107
 © Tetracore is a Registered Trademark



Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

Veterinary Real Time PCR Reagents Continued		
Product Name	Catalog Number	Product Description
<i>Bacillus anthracis</i> Lethal Factor (pXO1) PCR	TC-9020-YYY	Reagents for the specific detection of <i>Bacillus anthracis</i> Lethal Factor (PXO1)
<i>Bacillus anthracis</i> capA (pXO2) PCR	TC-9021-YYY	Reagents for the specific Detection of <i>Bacillus anthracis</i> capA (pXO2)
<i>Brucella spp.</i> PCR	TC-9022-YYY	Reagents for the specific detection of <i>Brucella spp.</i>
<i>Francisella tularensis</i> PCR	TC-9023-YYY	Reagents for the specific detection of <i>Francisella tularensis</i>
Japanese Encephalitis Virus RT-PCR	TC-9025-YYY	Reagents for the specific detection of Japanese Encephalitis Viral RNA
Orthopox Virus PCR	TC-9026-YYY	Reagents for the specific detection of Orthopox Viral DNA
Venezuelan Equine Encephalitis Virus RT-PCR	TC-9028-YYY	Reagents for the specific detection of Venezuelan Equine Encephalitis Viral RNA
Real Time PCR Primer/Probe Set (No Taq)		
Product Name	Catalog Number	Product Description
<i>Bacillus globigii</i> Primer/Probe Set	TC-8012-YYY	Primer/Probe set for the specific detection of <i>Bacillus globigii</i>
<i>Bacillus thuringiensis</i> subsp. <i>Kurstaki</i> Primer/Probe set	TC-8013-YYY	Primer/Probe set for the specific detection <i>Bacillus thuringiensis</i> subsp. <i>Kurstaki</i>
<i>Bacillus thuringiensis</i> subsp. <i>Israelensis</i> Primer/Probe set	TC-8014-YYY	Primer/Probe set for the specific detection of <i>Bacillus thuringiensis</i> subsp. <i>Israelensis</i>
Dengue Virus (Group; Serotypes 1,2,3,4) Primer/Probe set	TC-500X-YYY	Primer/Probe set for the specific detection of Dengue Viral RNA
Yellow Fever Virus Primer/Probe Set	TC-5006-YYY	Reagents for the specific detection of Yellow Fever Viral RNA

9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
 Phone: 240.268.5400 * Fax: 240.268.1107
 ® Tetracore is a Registered Trademark



Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

BioThreat Alert® Kits			
Product Name	Catalog Number	Product Description	Unit Size
Abrin BioThreat Alert® Kit	TC-8011-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (25 Strips)
BioThreat Alert® Sample Collection Kit	TC-8003-025	Materials used for obtaining and collecting samples to be used with BioThreat Alert Strips.	1 Re-Sealable Bag
Anthrax BioThreat Alert® Kit	TC-8004-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (25 Strips)
Bot Tox BioThreat Alert® Kit	TC-8005-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (25 Strips)
Brucella BioThreat Alert® Kit	TC-8002-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (25 Strips)
Orthopox BioThreat Alert® Kit	TC-8001-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (25 Strips)
Plague BioThreat Alert® Kit	TC-8006-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (25 Strips)
Proficiency BioThreat Alert® Kit	TC-8007-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (10 + Strips, 10 – Strips)
Ricin BioThreat Alert® Kit	TC-8008-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (25 Strips)
SEB BioThreat Alert® Kit	TC-8009-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (25 Strips)
Tularemia BioThreat Alert® Kit	TC-8010-025	Lateral Flow Immunoassay Test Strips: Includes Materials for Collection of 5 Samples	1 Box (25 Strips)

Rapid BioAlert Strip Reader (TC-3001-001)



9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
 Phone: 240.268.5400 * Fax: 240.268.1107
 ® Tetracore is a Registered Trademark



Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

BioThreat Alert® ELISA Kits				
Product Name	Catalog Number	Product Description	Test Format & Instrument Platforms	Unit Size
Abrin ELISA Kit	TC-4029-002 TC-4030-002 TC-4031-001 TC-4032-001	Sensitive Assays for the Specific Detection of Abrin in Complex Matrices	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/-) Pre-Coated (+/-)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate
<i>B. anthracis</i> ELISA Kit	TC-4005-002 TC-4006-002 TC-4007-001 TC-4008-001	Sensitive Assays for the Specific Detection of <i>B. anthracis</i> in Complex Matrices	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/only) Pre-Coated (+/only)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate
<i>B. globigii</i> ELISA Kit	TC-4034-002 TC-4036-002 TC-4033-001 TC-4035-001	Sensitive Assays for the Specific Detection of <i>B. globigii</i>	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/only) Pre-Coated (+/only)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate
<i>Botulinum</i> toxin A & B ELISA Kit	TC-4009-003 TC-4010-002 TC-4011-001 TC-4012-001	Sensitive Assays for the Specific Detection of Botulinum Toxin in Complex Matrices	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/only) Pre-Coated (+/only)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate
<i>Botulinum</i> Toxin A ELISA Kit	TC-4046-002 TC-4045-001	Sensitive Assays for the Specific Detection of <i>Botulinum</i> Toxin A	Non-Coated (+/-) Non-Coated (+/only)	Two 96 Well Plates One 96 Well Plate
<i>Botulinum</i> Toxin B ELISA Kit	TC-4039-002 TC-4037-002 TC-4040-001 TC-4038-001	Sensitive Assays for the Specific Detection of <i>Botulinum</i> Toxin B	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/only) Pre-Coated (+/only)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate
<i>Botulinum</i> Toxin C ELISA Kit	TC-4058-002 TC-4057-001	Sensitive Assays for the Specific Detection of <i>Botulinum</i> Toxin C	Non-Coated (+/-) Non-Coated (+/only)	Two 96 Well Plates One 96 Well Plate
<i>Botulinum</i> Toxin E ELISA Kit	TC-4054-002 TC-4053-001	Sensitive Assays for the Specific Detection of <i>Botulinum</i> Toxin E	Non-Coated (+/-) Non-Coated (+/only)	Two 96 Well Plates One 96 Well Plate
Brucella ELISA Kit	TC-4025-002 TC-4026-002 TC-4027-001 TC-4028-001	Sensitive Assays for the Specific Detection of Brucella In Complex Matrices	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/only) Pre-Coated (+/only)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate

9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
 Phone: 240.268.5400 * Fax: 240.268.1107
 ® Tetracore is a Registered Trademark



Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

BioThreat Alert® ELISA Kits				
Product Name	Catalog Number	Product Description	Test Format & Instrument Platforms	Unit Size
<i>F. tularensis</i> ELISA Kit	TC-4021-002 TC-4022-002 TC-4023-001 TC-4024-001	Sensitive Assays for the Specific Detection of <i>F. tularensis</i> in Complex Matrices	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/only) Pre-Coated (+/only)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate
MS2 Bacteriophage ELISA Kit	TC-4042-002 TC-4041-001	Sensitive Assays for the Specific Detection of MS2 bacteriophage	Non-Coated (+/-) Non-Coated (+/only)	Two 96 Well Plates One 96 Well Plate
Ricin ELISA Kit	TC-4001-002 TC-4002-002 TC-4003-001 TC-4004-001	Sensitive Assays for the Specific Detection of the Ricin	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/only) Pre-Coated (+/only)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate
SEB ELISA Kit	TC-4013-002 TC-4014-004 TC-4015-001 TC-4016-001	Sensitive Assays for the Specific Detection of SEB in Complex Matrices	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/only) Pre-Coated (+/only)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate
<i>Staphylococcal enterotoxin A</i> ELISA Kit	TC-4062-002 TC-4061-001	Sensitive Assays for the Specific Detection of <i>Staphylococcal enterotoxin A</i>	Non-Coated (+/-) Non-Coated (+/only)	Two 96 Well Plates One 96 Well Plate
<i>Y. Pestis</i> ELISA Kit	TC-4017-002 TC-4018-002 TC-4019-001 TC-4020-001	Sensitive Assays for the Specific Detection of <i>Y. pestis</i> in Complex Matrices	Non-Coated (+/-) Pre-Coated (+/-) Non-Coated (+/only) Pre-Coated (+/only)	Two 96 Well Plates Two 96 Well Plates One 96 Well Plate One 96 Well Plate

9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
 Phone: 240.268.5400 * Fax: 240.268.1107
 ® Tetracore is a Registered Trademark



Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

Antibody Products			
Product Name	Catalog Number	Product Description	Unit Size
Rabbit Anti-Abrin	TC-7027-001	Abrin Antibody	100 µg
	TC-7027-002		1 mg
	TC-7027-003		Other
Mouse Monoclonal Antibody to Abrin Isotype IgG21	TC-7026-001	Abrin Antibody	100 µg
	TC-7026-002		1 mg
	TC-7026-003		Other
Goat Anti- <i>B. anthracis</i> IgG	TC-7015-001	<i>B. anthracis</i> Antibody	100 µg
	TC-7015-002		1 mg
	TC-7015-003		Other
Rabbit Anti <i>B. anthracis</i> IgG	TC-7009-001	<i>B. anthracis</i> Antibody	100 µg
	TC-7009-002		1 mg
	TC-7009-003		Other
Goat Anti <i>B. globigii</i> IgG	TC-7014-01	<i>B. globigii</i> Antibody	100 µg
	TC-7014-02		1 mg
	TC-7014-03		Other
Rabbit Anti <i>B. globigii</i> IgG	TC-7008-001	<i>B. globigii</i> Antibody	100 µg
	TC-7008-002		1 mg
	TC-7008-003		Other
Mouse Monoclonal Antibody to <i>B. globigii</i> Isotype IgG	TC-7025-001	<i>B. globigii</i> Antibody	100 µg
	TC-7025-002		1 mg
	TC-7025-003		Other
Goat Anti-Botulinum Toxin A and B IgG	TC-7012-01	<i>Botulinum</i> Toxin A & B Antibody	100 µg
	TC-7012-02		1 mg
	TC-7012-03		Other
Rabbit Anti-Botulinum Toxin A and B IgG	TC-7007-001	<i>Botulinum</i> Toxin A & B Antibody	100 µg
	TC-7007-002		1 mg
	TC-7007-003		Other
Mouse Monoclonal Antibody to Botulinum Toxin A and B IgG	TC-7024-001	<i>Botulinum</i> Toxin A & B Antibody	100 µg
	TC-7024-002		1 mg
	TC-7024-003		Other
Mouse anti- <i>Botulinum</i> C Toxin, Isotype IgG1	TC-7037-001	<i>Botulinum</i> C Toxin Antibody	100 µg
	TC-7037-002		1 mg
	TC-7037-003		Other
Rabbit anti- <i>Botulinum</i> E Toxin	TC-7036-001	<i>Botulinum</i> E Toxin Antibody	100 µg
	TC-7036-002		1 mg
	TC-7036-003		Other
Rabbit anti- <i>Botulinum</i> F Toxin	TC-7035-001	<i>Botulinum</i> F Toxin Antibody	100 µg
	TC-7035-002		1 mg
	TC-7035-003		Other
Goat Anti-Brucella IgG1	TC-7011-01	Brucella Antibody	100 µg
	TC-7011-02		1 mg
	TC-7011-03		Other

9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
 Phone: 240.268.5400 * Fax: 240.268.1107
 © Tetracore is a Registered Trademark



Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

Antibody Products

Product Name	Catalog Number	Product Description	Unit Size
Mouse Monoclonal Antibody to Brucella Isotype IgG1	TC-7028-001 TC-7028-002 TC-7028-003	Brucella Antibody	100 µg 1 mg Other
Rabbit Anti- <i>E. herbicola</i> IgG	TC-7006-001 TC-7006-002 TC-7006-003	<i>E. herbicola</i> Antibody	100 µg 1 mg Other
Mouse Monoclonal Antibody to <i>E. herbicola</i> isotype IgG2b	TC-7029-001 TC-7029-002 TC-7029-003	<i>E. herbicola</i> Antibody	100 µg 1 mg Other
Rabbit Anti- <i>F. tularensis</i> IgG	TC-7005-001 TC-7005-02 TC-7005-03	<i>F. tularensis</i> Antibody	100 µg 1 mg Other
Mouse Monoclonal Antibody to <i>F. tularensis</i> Isotype IgG	TC-7021-001 TC-7021-002 TC-7021-003	<i>F. tularensis</i> Antibody	100 µg 1 mg Other
Rabbit Anti-MS2 IgG	TC-7004-001 TC-7004-002 TC-7004-003	MS2 Antibody	100 µg 1 mg Other
Mouse Monoclonal Antibody to MS2 Bacteriophage, Isotype IgG2a	TC-7033-001 TC-7033-002 TC-7033-003	MS2 Antibody	100 µg 1 mg Other
Rabbit Anti-Ovalbumin IgG	TC-7003-001 TC-7003-002 TC-7003-003	Ovalbumin Antibody	100 µg 1 mg Other
Mouse Monoclonal Antibody to Ovalbumin Isotype IgG1	TC-7034-001 TC-7034-002 TC-7034-003	Ovalbumin Antibody	100 µg 1 mg Other
Goat Anti-Ricin IgG	TC-7013-001 TC-7013-002 TC-7013-003	Ricin Antibody	100 µg 1 mg Other
Rabbit Anti Ricin IgG	TC-7002-001 TC-7002-002 TC-7002-03	Ricin Antibody	100 µg 1 mg Other
Mouse Monoclonal Antibody to Ricin Isotype IgG	TC-7017-001 TC-7017-002 TC-7017-002	Ricin Antibody	100 µg 1 mg Other
Goat Anti-SEB IgG	TC-7010-001 TC-7010-002 TC-7010-003	SEB Antibody	100 µg 1 mg Other
Rabbit Anti-SEB IgG	TC-7001-001 TC-7001-002 TC-7001-003	SEB Antibody	100 µg 1 mg Other

9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
 Phone: 240.268.5400 * Fax: 240.268.1107
 ® Tetracore is a Registered Trademark



Quick Reference Guide for Tetracore® Products
 Please visit our website at www.tetracore.com or contact
 Pete Pillay at (240) 268-5417 ppillay@tetracore.com

Antibody Products

Product Name	Catalog Number	Product Description	Unit Size
Mouse Monoclonal Antibody to SEB Isotype IgG ₁	TC-7020-001	SEB Antibody	100 µg
	TC-7020-002		1 mg
	TC-7020-003		Other
Mouse Monoclonal Antibody to SEC	TC-7031-001	SEC Antibody	100 µg
	TC-7031-002		1 mg
	TC-7031-003		Other
Rabbit Anti- <i>Y. pestis</i> IgG	TC-7022-001	<i>Y. pestis</i> Antibody	100 µg
	TC-7022-002		1 mg
	TC-7022-003		Other
Mouse Monoclonal Antibody to <i>Y. pestis</i> Isotype IgG	TC-7019-001	<i>Y. pestis</i> Antibody	100 µg
	TC-7019-002		1 mg
	TC-7019-003		Other

FDA Cleared Test for *B. Anthracis*

Clinical Products

Product Name	Catalog Number	Product Description	Test Format & Instrument Platforms	Unit Size
RedLine Alert™	TC-5123-010	Immunochromatographic test for the identification of <i>Bacillus</i> colonies cultured on sheep blood agar plates	NA	10 Cassettes/ 4 positive controls

Positive Result



Negative Result



9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
 Phone: 240.268.5400 * Fax: 240.268.1107
 © Tetracore is a Registered Trademark



Quick Reference Guide for Tetracore® Products
Please visit our website at www.tetracore.com or contact
Pete Pillay at (240) 268-5417 ppillay@tetracore.com

Instrument Products		
Product Name	Product Description	Catalog Number
Rapid BioAlert Strip Reader	A strip reader that can be used for TC-8XXX-XXX series products	TC-3001-001
Ultra Portable 4 Channel, Dual Color Real-Time PCR Device (Coming Soon)	A portable real-time PCR instrument that can be used for TC-9XXX-XXX series products	TC-3002-004

Ultra Portable 4-Well, Multiplex Capable (1 Color +IC or 2 Color) Real-Time PCR Device (TC-3002-004)



9901 Belward Campus Drive Ste 300, Rockville, Maryland 20850
Phone: 240.268.5400 * Fax: 240.268.1107
® Tetracore is a Registered Trademark